

DEVOLPE AMBU BAG



DEVICE DESCRIPTION

A bag valve mask, sometimes known by the proprietary name Ambu bag or generically as a manual resuscitator or "self-inflating bag", is a hand-held device commonly used to provide positive pressure ventilation to patients who are not breathing or not breathing adequately.

Features:

- ❖ Adult type 1500 ml
- ❖ Children type 450 ml
- ❖ Infant type 280 ml
- ❖ Consisting a set of Silicone Manual Resuscitator and resuscitation mask, 2 m PVC O2 safety tube and 2000 ml (for adult type) and 1600 ml PVC reservoir bags (for child & infant types).

TECHNICAL CHARACTERISTICS

Material: Plastic and its compositions

Ambu bag, reusable	Catalog code	Weight (kg)	Accessories
	AFAB/01/ABDP	275 gr	Available with safety valve (adult 60 mbar, child+infant 40mbar) and resuscitator masks for different patient types

Made in accordance with: ISO 13485:2016 standard

