



REF DUO-SLIDER STRETCHER BASE SUPPORT

Catalog number: SBS/02/DSDP

MANUFACTURER: DEVOLPE s.r.o.

REGISTERED TRADE MARK: DEVOLPE

REGISTERED ADDRESS: Nam. 14. října 1307/2, Prague 5, 150 00, Czech Republic

T/F: +420 605 700 078

Website: www.devolpe.com

Email: info@devolpe.com

Drafted in accordance with the standards:

EN 1865-1:2010-1



DEVICE DESCRIPTION

Duo-Slider Base Support is the device that is fitted to the floor and/or sidewall of an ambulance, on which a stretcher or incubator is located. The Device yields extensive solutions for paramedics to carry out Emergency Medical Services at the highest grade inside the ambulances together with its easiness of loading and unloading the stretcher on and off ambulances. It has sturdy frame, sliding options and space for spinal board or scoop stretcher. The ramp of the device is made by stainless steel and equipped with 2 gas springs-1 piece from each side to ensure the stability and smooth functionality of the ramp part.

New models are equipped with the blockage mechanism inside the tray for scoop stretcher/spinal board to prevent possible shakiness.

TECHNICAL CHARACTERISTICS

Material: Aluminum with stainless steel and steel compositions

Length (cm): 202 cm

Width (cm): 62 cm

Height (cm): 19 cm

Height of scoop stretcher/spinal board tray (cm): 9 cm

Slideability (cm): 20 cm

Weight (kg): 60 kg

Load capacity (kg): 250 kg



for more info please contact the manufacturer:

<https://www.devolpe.com/index.php?route=information/contact>



Devolve s.r.o.

Náměstí 14, Října 1307/2, Smíchov, 150 00 Prague 5, Czech Republic

+420 605 700 078 info@devolpe.com www.devolpe.com