

## DetaCarry Stretcher

**Catalogue code: MAS/02/DCDP**



### DEVICE DESCRIPTION

MonoCarry Stretcher (later called as “device”) is a device used to transport patients to hospitals for further medical treatment.

The device is operated by minimum of two users.

Key peculiarities of the device:

- It has 250 kg load capacity;
- Upper stretcher is detachable from the trolley by using locks;
- Gas springs equipped at both of backrest and footrest parts ease the positioning and provides unlimited positions during use;
- The device has four 360° swivel wheels which two of them equipped with brakes;
- Height adjustment system has safety mechanism which prevents any adverse event during operations and positioning and secures the safety of users and patients;

Features:

- The device paint is water, disinfectants and scratch resistant;
- The frame consists of sturdy lightweight non twisting construction enabling use of cardiopulmonary resuscitation;
- The backrest can be raised the at least 75° and more than that;
- The device is equipped with the sets of 3 quick-release patient restraints;



- Lying part of the device is produced by antibacterial, antifungal, anti-stain and anti-pudrid that easing to wash, clean all kinds of liquids and to disinfect;
- The frame is made non-twisted and the device is free of edges;
- The mattress of the device is produced with anti-bacterial material which is easy to clean and disinfect.

Used materials: Aluminum, stainless steel, steel and their composites

## TECHNICAL CHARACTERISTICS

Product	Catalogue code	Length (cm)	Width (cm)	Height (cm)	Dimensions of the main carriage wheels	Maximum load height (cm)	Weight (kg)	Load capacity (kg)	Positions
MonoCarry stretcher	MAS/01/MCDP	195	58	102	200 mm	70 cm (can be customized)	45	250	Trendelenburg; Reverse-Trendelenburg, Middle position; Load position; Floor position

Made in accordance with EN 1865-1 standard

And under certified ISO 13485:2016 standard

## Accessories

- ❖ intravenous pole (IV pole);
- ❖ quick-release system cross model (can be chosen or straight model is provided automatically);
- ❖ stretcher fixation system;
- ❖ Oxygen cylinder holder;
- ❖ Mattress (as extra spare part)

**\*note:**

**for more info please contact to the manufacturer. Photos are followings:**

