

EXTRICATION COLLAR



REF

Catalog number: ID/02/ECDP

Manufacturer: Devolpe s.r.o.

DEVICE DESCRIPTION

is a medical device used to support a person's neck. It is also applied by emergency personnel to those who have had traumatic head or neck injuries, and can be used to treat chronic medical conditions. They can also be used for extricating patients from car accidents who have suspected whiplash, in such cases the device is used until X-rays determine any fracture. The device can be used with the combination of Kendrick Extrication Device or backboard. Additionally, X-ray, CT and MRI procedures can be performed with the device in place.

TECHNICAL CHARACTERISTICS

Product	Material	Catalog code	Length (cm)	Width (cm)	Height (cm)	Weight (kg)	Size range
Extrication Collar	Hard polyethylene material	ID/02/ECDP	54	17	1.3	0.15	Universal combining 4 sizes, S, M, L and XL

Made in accordance with: ISO 13485:2016 standard

for more info please contact the manufacturer:

<https://www.devolpe.com/index.php?route=information/contact>

