

Escachair with stair belts



REF

Catalog number: FPTC/02/ECDP

Manufacturer: Devolpe s.r.o.

DEVICE DESCRIPTION

Escachair with stair belts (later called as “device”) is a device used to transport patients to hospitals for further medical treatment in a sitting position.

The device is operated by minimum of two users.

Key peculiarities of the device:

- It has 220 kg load capacity;
- Sitting part can be made of a metal or aluminum alloy covered by fabric, depending on the order request;
- The device has four 360° swivel wheels which two back wheels equipped with brakes;
- The device has pair of non-slip handles at the foot bar and top frame of the device;
- Sitting part and backrest is equipped with the foam pads connected with the high quality velcro to the device to ensure the comfortable transport of the patient;
- Footrest provides comfortness during transportation.

Features:

- The device paint is water, desinfectants and scratch resistant;
- The frame consists of sturdy lightweight non twisting construction that is free of edges;



Devolpe s.r.o.

📍 Náměstí 14. Října 1307/2, Smíchov, 150 00 Prague 5, Czech Republic
📞 +420 605 700 078 ✉ info@devolpe.com 🌐 www.devolpe.com

- The device is equipped with the sets of 2 quick-release patient restraints;
- The pads covering sitting part and backrest of the device are produced by antibacterial, antifungal, anti-stain and anti-pudrid material that is waterproof, easy to wash and clean all kinds of liquids & petrol-oil and to desinfect.

Used materials: Aluminum, stainless steel, steel and their composites and fabric.

TECHNICAL CHARACTERISTICS

Product	Catalogue code	Length (cm)	Width (cm)	Height (cm)	Dimensions of the rear wheels (mm)	Weight (kg)	Load capacity (kg)
EscaChair with stair belts	FPTC/02/ECDP	75 (non-folded) 30 (folded)	48	105	100	*10 kg (without optional parts)	220

**Note: Stair Belts and their parts of the chair are the optional parts of the device, hence, the weight of the device with the optional parts is 19 kg.*

Made in accordance with EN 1865-1 standard

And under certified ISO 13485:2016 standard

Accessories:



- Quick-release patient restraints
Catalog code: FPTCA/01/DP



- Sitting & Lying pads
Catalog code: FPTCA/02/DP

for more info please contact the manufacturer:

<https://www.devolpe.com/index.php?route=information/contact>



Devolve s.r.o.

Náměstí 14. Října 1307/2, Smíchov, 150 00 Prague 5, Czech Republic
+420 605 700 078 info@devolve.com www.devolpe.com

