

FOLDABLE STRETCHER IN 2 FOLDS



Product`s description and intended usage:

This product is designated transferring patients at pre-hospital immobilizations under hard-accessible conditions.

It is made of high-quality aluminum and retardant fabric oxford for lying part. The device is folded into 2 folds before and after usage. The product has restraint system for patient`s safety.

As the sole manufacturer, we can change the color, material of the lying part and some specifications if requested.

TECHNICAL SPECIFICATIONS:

Product	Catalogue code	Unfolded Length (cm)	Unfolded Width (cm)	Unfolded Height (cm)	Folded Length (cm)	Folded Width (cm)	Folded Height (cm)	Weight (kg)	Load capacity (kg)
Foldable stretcher into 2 Folds	PUS/01/FSDP	215	58	12	108	26	17	6	200

Made in accordance with: ISO 13485:2016 standard &

EN 1865-1 standard

Class category is Class I under Medical Device Regulation of EU.



DEVOLPE
MEDICAL

Accessories

The product is available with its carrying bag.

