



DEVOLPE
MEDICAL

FOLDABLE STRETCHER IN 4 FOLDS CIVIL & MILITARY TYPE



Product`s description and intended usage:

This product is designated transferring patients at pre-hospital immobilizations under hard-accessible conditions.

It is made of high-quality aluminum and retardant fabric oxford for lying part. The device is folded into 4 folds before and after usage. The product has 2 pcs. of restraint system for patient`s safety.

As the sole manufacturer, we can change the color, material of the lying part and some specifications if



requested.

TECHNICAL SPECIFICATIONS:

Product	Catalogue code	Unfolded Length (cm)	Unfolded Width (cm)	Unfolded Height (cm)	Folded Length (cm)	Folded Width (cm)	Folded Height (cm)	Weight (kg)	Load capacity (kg)
Foldable stretcher into 4 Folds	PUS/01/FSDP	215	55	15	50	25	18	9	250

Made in accordance with: ISO 13485:2016 standard &

EN 1865-1 standard

Class category is Class I under Medical Device Regulation of EU.

Accessories

The product is available with its carrying bag. Dimensions are 56 x 28 x 20 cm



Photo from another military color

