



## KENDRICK EXTRICATION DEVICE

# REF

Catalog number: ID/01/KEDDP

MANUFACTURER: DEVOLPE s.r.o.

REGISTERED TRADE MARK: DEVOLPE

REGISTERED ADDRESS: Nam. 14. října 1307/2, Prague 5, 150 00, Czech Republic

T/F: +420 605 700 078

Website: [www.devolpe.com](http://www.devolpe.com)

Email: [info@devolpe.com](mailto:info@devolpe.com)

*Drafted in accordance with the standards:*

*EN 1865-1:2010-1*



Devolpe s.r.o.

Náměstí 14. Října 1307/2, Smíchov, 150 00 Prague 5, Czech Republic  
+420 605 700 078 [info@devolpe.com](mailto:info@devolpe.com) [www.devolpe.com](http://www.devolpe.com)

## DEVICE DESCRIPTION

As the kind of immobilization device, Kendrick Extrication devices are used for extricating and immobilizing patients from motor vehicles in case of accidents. They can also be used for extricating patients from other confined spaces safely. The semi-rigidity of the device secures the head and cervico-thoracic area in a neutral position which minimizes the risk of further injury. Additionally, X-ray, CT and MRI procedures can be performed with the device in place.

## TECHNICAL CHARACTERISTICS

Material: Aluminum with stainless steel and steel compositions

Product	Material	Catalog code	Length (cm)	Width (cm)	Height (cm)	Weight (kg)
Kendrick Extrication Device (KED)	PVC material and composites	ID/01/K EDDP	88	76	3	2

Made in accordance with: EN 1865-1 standard

## Optional Accessories



## Carrying bag for KED

for more info please contact the manufacturer:

<https://www.devolpe.com/index.php?route=information/contact>



Devolve s.r.o.

Náměstí 14. Října 1307/2, Smíchov, 150 00 Prague 5, Czech Republic

+420 605 700 078 [info@devolve.com](mailto:info@devolve.com) [www.devolpe.com](http://www.devolpe.com)