

## MANUAL SUCTION PUMP



### DEVICE DESCRIPTION

Manual Suction Pump is a light and portable suction unit that can be operated by one hand. Manual Hand Suction Device features a disposable canister that holds aspirator contents securely. Portable Manual Suction Pump is designed to provide simple operation and maintenance in all emergency suction applications and is perfect for emergency suction of airways secretions with safety, quickness, and efficacy. Each unit includes a reusable pump, a suction canister and tubing. Portable Suction Pump is designed with built-in overflow protection.

Features:

- ❖ High vacuum, light , portable, emergency suction device
- ❖ Built in 3M HEPA filter for bacterial and viral protection
- ❖ Adult and child suction settings
- ❖ Built in overflow protection seal
- ❖ 100°swivel allows better accessibility in tight situations
- ❖ Efficient operation and simple maintenance
- ❖ One hand operation
- ❖ Large volume capacity

### TECHNICAL CHARACTERISTICS

Material: Plastic and its compositions

Made in accordance with: ISO 13485:2016 standard

Manual Suction pump, reusable	Catalog code	Length (cm)	Height (cm)	Weight (kg)	Included
	AFAB/02/MP DP	18	7	275 gr	1 reusable handheld suction pump unit. 1 disposable collection container (300mL). 1 suction catheter for adult (OD 17mm). 1 suction catheter for child (OD .8mm).

