



DEVOLPE
MEDICAL





Devolpe s.r.o
Nam. 14. Rijna 1307/2, 150 00, Prague 5
ICO: 09776885, DICO: CZ09776885

| | | |
|---|--|---|
| <p>1-piece Medical pouch including accessories below</p> <p>Manufactured by Devolpe Medikal Saglik Urunleri Ic ve Dis (Turkish division of Devolpe s.r.o.)</p> | <p>Made of high durable 500D Cordura (nylon also can be chosen according to the buyer`s needs) dimensions-20 cm length, 15 cm width, 10 cm height, weight 0.350 kg</p> |  |
| <p>1-piece Sterile gauze bandage</p> <p>Manufactured by Devolpe Medikal Saglik Urunleri Ic ve Dis (Turkish division of Devolpe s.r.o.)</p> | <p>gauze type 20 (density 27 g/m2), Clarity within length and width variations, ease of unwinding (the threads of the edges are not tangled). High capillarity (about 11 cm), breaking load on the length above the norm (about 90 N), high degree of whiteness, dimensions-7 m length, 14 cm width, with printed instructions in Ukrainian language</p> |  |
| <p>1-piece Emergency bandage/Israeli bandage</p> <p>Manufactured by Devolpe Medikal Saglik Urunleri Ic ve Dis (Turkish division of Devolpe s.r.o.)</p> | <p>The Emergency Bandage consists of a sterile, non-adherent pad that applies pressure to any site and elastic strap, can be easily wrapped and secured, and has an additional application, like a tourniquet, to further constrict blood- flow. Dimensions are 15cm x 7cm (4inches) with printed instructions in Ukrainian language</p> |  |
| <p>1-piece Adhesive plaster</p> <p>Manufactured by Devolpe Medikal Saglik Urunleri Ic ve Dis (Turkish division of Devolpe s.r.o.)</p> | <p>Adhesive Strapping Tape-Dimensions 2.5 cm x 5 m (pack of 5)</p> |  |
| <p>1-piece Isothermal blanket</p> <p>Manufactured by Devolpe Medikal Saglik Urunleri Ic ve Dis (Turkish division of Devolpe s.r.o.)</p> | <p>Material: Plastic, metallized with aluminum – Measures: 200 x 210 cm, Color: Gold/Silver, with printed instructions in Ukrainian language</p> |  |



DEVOLPE
MEDICAL

Devolpe s.r.o
Nam. 14. Rijna 1307/2, 150 00, Prague 5
ICO: 09776885, DICO: CZ09776885

| | | |
|---|--|---|
| <p>2-piece Tactical tourniquet</p> <p>Manufactured by Devolpe Medikal Saglik Urunleri Ic ve Dis (Turkish division of Devolpe s.r.o.)</p> | <p>Made of high quality plastic with Velcro, dimensions – 95 cm length, 3.8 cm widths, weight 90 gr, with printed instructions in Ukrainian language</p> |  |
| <p>1-piece Scissors</p> <p>Manufactured by Devolpe Medikal Saglik Urunleri Ic ve Dis (Turkish division of Devolpe s.r.o.)</p> | <p>Designed for cutting the clothes and shoes, material is stainless steel, 19 cm length (7,5 inches), with printed instructions in Ukrainian language on the package</p> |  |
| <p>1-piece Nasopharyngeal airway with lubricants</p> <p>Manufactured by Devolpe Medikal Saglik Urunleri Ic ve Dis (Turkish division of Devolpe s.r.o.)</p> | <p>Nasopharyngeal airway, material PVC, packed in separate packs, sizes: 2.5mm-9.0mm, length is 170 mm, with printed instructions in Ukrainian language</p> |  |
| <p>1-piece Blue marker medical (surgical)</p> <p>Manufactured by Devolpe Medikal Saglik Urunleri Ic ve Dis (Turkish division of Devolpe s.r.o.)</p> | <p>Leather marker with round tip, contains permanent ink with low odor, track width: 1mm with ventilation cap</p> |  |
| <p>3-pair Nitrile gloves</p> <p>Manufactured by Devolpe Medikal Saglik Urunleri Ic ve Dis (Turkish division of Devolpe s.r.o.)</p> | <p>Blue nitrile gloves, medical, powder-free and latex-free, size XL proportionally divided in each pouch and packed in paper + film package with printed instructions in Ukrainian language (3 pairs in each pair packed in separate package)</p> |  |
| <p>1-piece Hemostatic tamponing bandage Medclot brand</p> | <p>It has Hemostatic agent, 7.5 x 370 cm dimensions</p> |  |
| <p>1-piece The Fox eye shield</p> <p>Manufactured by Devolpe Medikal Saglik Urunleri Ic ve Dis (Turkish division of Devolpe s.r.o.)</p> | <p>A protective pad for an injured or operated eye. It is fixed thanks to the adhesive base. Dimensions: 7.3 cm x 6 cm Material: aluminum, fabric</p> |  |