

MonoCarry Stretcher

Catalogue code: MAS/01/MCDP



DEVICE DESCRIPTION

MonoCarry Stretcher (later called as “device”) is a device used to transport patients to hospitals for further medical treatment.

The device is operated by minimum of two users.

Key peculiarities of the device:

- It has 250 kg load capacity;
- Gas springs equipped at both of backrest and footrest parts ease the positioning and provides unlimited positions during use;
- The device has four 360° swivel wheels which two of them equipped with brakes;
- Height adjustment system has safety mechanism which prevents any adverse event during operations and positioning and secures the safety of users and patients;

Features:

- The device paint is water, disinfectants and scratch resistant;
- The frame consists of sturdy lightweight non twisting construction enabling use of cardiopulmonary resuscitation;
- The backrest can be raised the at least 75° and more than that with the support of gas spring;
- The legrest can be adjusted with the support of gas spring too;
- The device is equipped with swing-down side rails;
- The device is equipped with the sets of 3 quick-release patient restraints;



DEVOLPE
MEDICAL

- Lying part of the device is produced by antibacterial, antifungal, anti-stain and anti-pudrid that easing to wash, clean all kinds of liquids and to disinfect;
- The frame is made non-twisted and the device is free of edges;
- The mattress of the device is produced with anti-bacterial material which is easy to clean and disinfect.

Used materials: Aluminum, stainless steel, steel and their composites

TECHNICAL CHARACTERISTICS

Product	Catalogue code	Length (cm)	Width (cm)	Height (cm)	Dimensions of the main carriage wheels	Maximum load height (cm)	Weight (kg)	Load capacity (kg)	Positions
MonoCarry stretcher	MAS/01/MCDP	196	56	88	200 mm	70 cm (can be customized)	45 kg 53 kg (with accessories)	250	Trendelenburg; Reverse-Trendelenburg, Middle position; Load position; Floor position

Made in accordance with EN 1865-1 standard

And under certified ISO 13485:2016 standard

Accessories

- ❖ intravenous pole (IV pole);
- ❖ quick-release system cross model (can be chosen or straight model is provided automatically);
- ❖ stretcher fixation system;
- ❖ Oxygen cylinder holder;
- ❖ Mattress (as extra spare part)


***note:**

for more info please contact to the manufacturer. Photos are followings:





DEVOLPE
MEDICAL





DEVOLPE
MEDICAL

