



**REF**

**SCOOP STRETCHER**

**Catalog number: PUS/SSDP**

**MANUFACTURER: DEVOLPE s.r.o.**

**REGISTERED TRADE MARK: DEVOLPE**

**REGISTERED ADDRESS: Nam. 14. října 1307/2, Prague 5, 150 00, Czech Republic**

**T/F: +420 605 700 078**

**Website: [www.devolpe.com](http://www.devolpe.com)**

**Email: [info@devolpe.com](mailto:info@devolpe.com)**

***Drafted in accordance with the standards:***

***EN 1865-1:2010-1***



## DEVICE DESCRIPTION

Scoop stretcher is a device for displacement from a ground level to fixed stretcher of the injured suspicious of limb fractures or displacement of patients from a surface to the other.

Material: Aluminum alloy and composites

## TECHNICAL CHARACTERISTICS

Unfolded length (cm)	202
Unfolded width (cm)	42
Unfolded height (cm)	7
Folded length (cm)	118
Folded width (cm)	42
Folded height (cm)	7
Weight (kg)	8 kg
Load capacity (kg)	250

## Accessories



### Spider straps

Catalog code: ASB/01/SSDP

Material: heavy grade polypropylene or nylon webbing

### Carrying bag

Catalog code: ASB/02/CBDP

Material: polyester 600D



for more info please contact the manufacturer:

<https://www.devolpe.com/index.php?route=information/contact>



Devolve s.r.o.

Náměstí 14. Října 1307/2, Smíchov, 150 00 Prague 5, Czech Republic  
+420 605 700 078 info@devolpe.com www.devolpe.com