



DEVOLPE
MEDICAL

SPINE BOARD



REF

Catalog number: LS/01/SBDP

MANUFACTURER: DEVOLPE s.r.o.

REGISTERED TRADE MARK: DEVOLPE

REGISTERED ADDRESS: Nam. 14. října 1307/2, Prague 5, 150 00, Czech Republic

T/F: +420 605 700 078

Website: www.devolpe.com

Email: info@devolpe.com

Drafted in accordance with the standards:

EN 1865-1:2010-1



Devolpe s.r.o.

Náměstí 14. Října 1307/2, Smíchov, 150 00 Prague 5, Czech Republic

+420 605 700 078 info@devolpe.com www.devolpe.com



DEVICE DESCRIPTION

This product is designated transferring and transporting the patients at pre-hospital immobilizations and conditions. It is made high quality HDPE material and equipped with 8 pins. It has floating capacity too.


TECHNICAL CHARACTERISTICS

Material: Aluminum with stainless steel and steel compositions

| Product | Type | Catalog code | Length (cm) | Width (cm) | Height (cm) | Weight (kg) | Load capacity (kg) |
|-------------|-------|--------------|-------------|------------|-------------|-------------|--------------------|
| Spine Board | Adult | LS/01/SBDP | 184 | 44 | 5.5 | 5.5 | 250 kg |

Made in accordance with: EN 1865-1 standard

Accessories



Spider straps

Catalog code: ASB/01/SSDP

Material: heavy grade polypropylene or nylon webbing

