



**DEVOLPE**  
MEDICAL

## VACUUM MATTRESS



## REF

**Catalog number: VM/02/VMPDP**

**MANUFACTURER: DEVOLPE s.r.o.**

**REGISTERED TRADE MARK: DEVOLPE**

**REGISTERED ADDRESS: Nam. 14. října 1307/2, Prague 5, 150 00, Czech Republic**

**T/F: +420 605 700 078**

**Website: [www.devolpe.com](http://www.devolpe.com)**

**Email: [info@devolpe.com](mailto:info@devolpe.com)**

***Drafted in accordance with the standards:***

***EN 1865-1:2010-1***



Devolpe s.r.o.

Náměstí 14. Října 1307/2, Smíchov, 150 00 Prague 5, Czech Republic

+420 605 700 078 [info@devolpe.com](mailto:info@devolpe.com) [www.devolpe.com](http://www.devolpe.com)



## DEVICE DESCRIPTION

Vacuum mattress is a full spine immobilization device for multi-trauma patients with suspected spinal injuries. It is highly useful when the patient needs immobilization for more than 1 hour. The vacuum system allows soft and rigid immobilization, both useful on emergency transport vehicles.

Material: material combination with PU, PVC, EVA and new plastic materials.

## TECHNICAL CHARACTERISTICS

Product	Type	Models	Catalog code	Length (cm)	Width (cm)	Height (cm)	Weight (kg)
Vacuum Mattress	Adult	Solo-chamber	VM/01/VMSDP	200	80	3	7.95
Vacuum Mattress	Adult	Poli-chamber	VM/02/VMPDP	200	80	3	7.95
Vacuum Mattress (longer version)	Adult	Solo-chamber	VM/01/VMLSDP	210	80	3	7.95
Vacuum Mattress (longer version)	Adult	Poli-chamber	VM/02/VMPLDP	210	80	3	7.95

Made in accordance with: ISO 13485:2016 standard &

EN 1865-1 standard

## Accessories



**Manual Pump**  
**For vacuum products**



**Vacuum mattress`**  
**bag**

